



Vitrex Serum Clot Activator Tubes

Vitrex Serum Clot Activator Tubes are used for clinical chemistry, immunology, protein electrophoresis, serology, microbiology and toxicology. The inner wall of the Serum Clot Activator tube is specially coated with clot activator, which activates the coagulation process.



Transit packaging

Pallet

Packaging:

Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.
Serum CAT/Clot Activator	1	13 x 75	Serum Clot Activator	Paper	101201+ *
Serum CAT/Clot Activator	2	13 x 75	Serum Clot Activator	Paper	102202+ *
Serum CAT/Clot Activator	3	13 x 75	Serum Clot Activator	Paper	103203+ *
Serum CAT/Clot Activator	3	13 x 100	Serum Clot Activator	Paper	108203+ **
Serum CAT/Clot Activator	4	13 x 75	Serum Clot Activator	Paper	104204+ *
Serum CAT/Clot Activator	4	13 x 100	Serum Clot Activator	Paper	108204+ **
Serum CAT/Clot Activator	5	13 x 75	Serum Clot Activator	Paper	109205+ *
Serum CAT/Clot Activator	5	13 x 100	Serum Clot Activator	Paper	108205+ **
Serum CAT/Clot Activator	6	13 x 100	Serum Clot Activator	Paper	105206+ **
Serum CAT/Clot Activator	7	13 x 100	Serum Clot Activator	Paper	108207+ **
Serum CAT/Clot Activator	8	16 x 100	Serum Clot Activator	Paper	108208+ ***
Serum CAT/Clot Activator	9	16 x 100	Serum Clot Activator	Paper	106209+ ***
Serum CAT/Clot Activator	10	16 x 100	Serum Clot Activator	Paper	107200+ ***

Shelf life 18 months