



Vitrex EDTA Tubes

Vitrex EDTA Tubes are used for whole blood analysis within haematology. The range of application is: hematology, PCR analysis, blood typing. EDTA and Gel is mainly used for molecular diagnostics. The gel forms a resistant barrier between the cellular components of blood and serum during centrifugation.

K2EDTA

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Plain	+ Gel
					Article No.	Article No.
K2E/K2EDTA	1	13 x 75	K2EDTA	Paper	101901+ *	112201+ *
K2E/K2EDTA	2	13 x 75	K2EDTA	Paper	102902+ *	112202+ *
K2E/K2EDTA	3	13 x 75	K2EDTA	Paper	103903+ *	111203+ *
K2E/K2EDTA	3	13 x 100	K2EDTA	Paper	103930+ **	113203+ **
K2E/K2EDTA	3.5	13 x 75	K2EDTA	Paper		112235+ *
K2E/K2EDTA	4	13 x 75	K2EDTA	Paper	104904+ *	112240+ *
K2E/K2EDTA	4	13 x 100	K2EDTA	Paper	105904+ **	112204+ **
K2E/K2EDTA	5	13 x 75	K2EDTA	Paper	105905+ *	
K2E/K2EDTA	5	13 x 100	K2EDTA	Paper	105950+ **	112205+ **
K2E/K2EDTA	6	13 x 100	K2EDTA	Paper	106906+ **	112206+ **
K2E/K2EDTA	7	13 x 100	K2EDTA	Paper	106907+ **	
K2E/K2EDTA	8	16 x 100	K2EDTA	Paper	106908+ ***	113208+ ***
K2E/K2EDTA	8.5	16 x 100	K2EDTA	Paper		113209+ ***
K2E/K2EDTA	9	16 x 100	K2EDTA	Paper	107909+ ***	113290+ ***
K2E/K2EDTA	9.5	16 x 100	K2EDTA	Paper	108909+ ***	

Shelf life 18 months



Transit packaging

Pallet

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

K3EDTA

Plain + Gel

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.	Article No.
K3E/K3EDTA	1	13 x 75	K3EDTA	Paper	111101+ *	122101+ *
K3E/K3EDTA	2	13 x 75	K3EDTA	Paper	112102+ *	122102+ *
K3E/K3EDTA	3	13 x 75	K3EDTA	Paper	113103+ *	123103+ *
K3E/K3EDTA	3	13 x 100	K3EDTA	Paper	113130+ **	123130+ **
K3E/K3EDTA	3.5	13 x 75	K3EDTA	Paper		124135+ *
K3E/K3EDTA	4	13 x 75	K3EDTA	Paper	113104+ *	125140+ *
K3E/K3EDTA	4	13 x 100	K3EDTA	Paper	114104+ **	125104+ **
K3E/K3EDTA	5	13 x 75	K3EDTA	Paper	114105+ *	
K3E/K3EDTA	5	13 x 100	K3EDTA	Paper	114150+ **	125150+ **
K3E/K3EDTA	6	13 x 100	K3EDTA	Paper	116106+ **	126106+ **
K3E/K3EDTA	7	13 x 100	K3EDTA	Paper	116107+ **	
K3E/K3EDTA	8	16 x 100	K3EDTA	Paper	114108+ ***	127108+ ***
K3E/K3EDTA	8.5	16 x 100	K3EDTA	Paper	114185+ ***	128185+ ***
K3E/K3EDTA	9	16 x 100	K3EDTA	Paper	117109+ ***	129109+ ***

Shelf life 18 months