



## Chiraflex® IV Cannulas and Chiraflex® Safety IV Cannulas

Chiraflex IV Cannulas and Chiraflex Safety IV Cannulas are intended for peripheral vascular application. The catheter is made from Teflon material and available in Safety. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. They are packed in boxes of 50 pcs, cartons of 1000 pcs.

**USER CONVENIENCE:** Compatible with all standard medical devices for wide application.

**USER SAFETY:** Safety cannula is available.

**PATIENT SAFETY:** Single-use only. Sterilised by ETO.

Colour Code	Size G/™	Size Metric	Article No. With Port With Wings	Article No. Without Port With Wings	Article No. Without Port Without Wings	Article No. SAFETY With Port With Wings
Purple	26G	0.60 x 19	CHFX0102619	CHFX0202619	CHFX0302619	CHFXS0102619
Yellow	24G	0.70 x 19	CHFX0102419	CHFX0202419	CHFX0302419	CHFXS0102419
Blue	22G	0.90 x 25	CHFX0102225	CHFX0202225	CHFX0302225	CHFXS0102225
Pink	20G	1.10 x 32	CHFX0102032	CHFX0202032	CHFX0302032	CHFXS0102032
Green	18G	1.30 x 45	CHFX0101845	CHFX0201845	CHFX0301845	CHFXS0101845
White	17G	1.50 x 45	CHFX0101745	CHFX0201745	CHFX0301745	CHFXS0101745
Grey	16G	1.70 x 45	CHFX0101645	CHFX0201645	CHFX0301645	CHFXS0101645
Orange	14G	2.10 x 45	CHFX0101445	CHFX0201445	CHFX0301445	CHFXS0101445

Shelf life 5 years