



## Chiraplus™ Infusion set

Chiraplus Infusion sets include 15 µm filter, ventilation port and Luer lock adapter. Chiraplus P is intended for pressure infusion and does not include flash ball. Chiraplus G is intended for gravity infusion and includes flash ball.

Product	Chamber Vol. (ml)	Tube length (mm)	Packaging Type	Pcs./Carton	Article No.
Chiraplus P	10.5	1500	Blister	100	CHP336
Chiraplus P Dehp free	10.5	1500	Blister	100	CHP536
Chiraplus P 12 Dehp free	12	1500	Blister	100	CHP506
Chiraplus P 190	10.5	1900	Blister	100	CHP339
Chiraplus P 190 Safety	10.5	1900	Blister	100	CHP338
Chiraplus P Dehp free Autofusion 180	15	1800	Blister	300	CHP576
Chiraplus G	10.5	1500	Blister	100	CHG331
Chiraplus G	10.5	1500	Poly bag	500	CHG317
Chiraplus G Dehp free	10.5	1500	Blister	100	CHG531
Chiraplus G 12	12	1500	Blister	100	CHG301
Chiraplus G 12 Dehp free	12	1500	Blister	100	CHG501
Chiraplus G	8.5	1500	Poly bag	500	CHG318
Chiraplus G Dehp free	8.5	1500	Poly bag	500	CHG518
Chiraplus G with needle	8.5	1500	Poly bag	500	CHG319
Chiraplus G with metal spike	8.5	1500	Poly bag	500	CHG316
Chiraplus G PVC free	12	1500	Blister	100	CHG731
Chiraplus G PVC free opaque	12	1900	Blister	100	CHG901

Shelf life 5 years