



Vitrex Blood Gas Tubes

Vitrex capillary tubes for blood gas analyses are made from heavy walled glass, suitable for 2.3 mm Ø adaptors. The uniform heparinisation secure proper function of the tubes for blood gas analysis. Vitrex Blood Gas Tubes are also available in plastic with sodium heparin in the same sizes as Vitrex Blood Gas Tubes in glass.

Accessories

- Stirrers packed in vials of 250 pcs.
- End Caps packed in boxes of 500 pcs.
- Magnets packed in bags of 10 pcs.





Blood Gas Tubes

| Vol. (µl) | Pcs. | d (mm) +/- 0.05 mm | D (mm) +/- 0.05 mm | L (mm) +/- 1 mm | Na Hep. Article No. | Li Hep. Article No. |
|-----------|---------|--------------------|--------------------|-----------------|---------------------|---------------------|
| 130 (min) | 4 x 250 | 1.60 | 2.30 | 75 | 182513 | 182613 |
| 175 (min) | 4 x 250 | 1.60 | 2.30 | 100 | 182413 | 183313 |
| 220 (min) | 4 x 250 | 1.60 | 2.30 | 125 | 182313 | 182013 |

Special types upon request.

Shelf life 4 years

Plastic Blood Gas Tubes

| Vol. (µl) | Pcs. | d (mm) +/- 0.05 mm | D (mm) +/- 0.05 mm | L (mm) +/- 1 mm | Plastic Na Hep. Article No. |
|-----------|---------|--------------------|--------------------|-----------------|-----------------------------|
| 130 (min) | 4 x 200 | 1.60 | 2.30 | 75 | 282513 |
| 175 (min) | 4 x 200 | 1.60 | 2.30 | 100 | 282413 |
| 220 (min) | 4 x 200 | 1.60 | 2.30 | 125 | 282313 |

Other types available upon request.

Shelf life 2 years

Accessories

| Description | For tubes with D | Pcs. | Article No. |
|---------------------|------------------|----------|-------------|
| End caps for D (mm) | 2.20 - 2.80 mm | 10 x 500 | 140813 |
| End caps for D (mm) | 1.50 - 2.20 mm | 10 x 500 | 140713 |
| Stirrers | | 10 x 250 | 140913 |
| Magnets | | 1 x 10 | 142400 |

