



Vitrex Flex 3

Depth Adjustable Safety lancet

Vitrex Flex 3 safety lancet is an All-in-One lancet for emergency use due to its adjustable design, which includes three different penetration depths for a single unit. The lancet is intended for blood glucose and other in vitro diagnostic tests.

USER CONVENIENCE: Pre-loaded and pressure activated design. The trigger position indicates the before and after use.

USER SAFETY: Needle is fully shielded before and after use.

COMFORT: High quality tri-bevel needle with high speed penetration.

PATIENT SAFETY: Single-use only. Sterilised by gamma irradiation.



Avoid misfiring when adjusting the depth. Ergonomic design for improved holding.



Improved safety cap removal – cap automatically semi ejects, when turned to remove.

Technical Specifications

Colour Code	Specification	Penetration Depth	Article No.
Blue 	23G	1.30 mm - 1.80 mm - 2.30 mm	395813
Yellow 	26G	1.30 mm - 1.80 mm - 2.30 mm	396813
Orange 	28G	1.30 mm - 1.80 mm - 2.30 mm	397813

Shelf life 5 years




Transit packaging

Pallet

Packaging:

Primary packaging: 100 pcs.
 Transit packaging: 1,000 pcs.
 Pallet: 84,000 pcs.



Penetration depth adjustment is done without touching the activation mechanism.