



Chiravac™ Blood Collection needles

Chiravac Blood Collection needles are made of polypropylene hub and stainless steel tube. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The needles are intended for blood collection with vacuum tubes. They are individually packed in boxes of 100 pcs, cartons of 1000 pcs. Upon request the needles are also available as safety or invisible version.

USER CONVENIENCE: Colour-coded needle hub, screw-on needle hub compatible with all available holders, easy-to-open packaging and thin wall tube material.

COMFORT: Grinded, precision lancet point and silicone treated tip for smooth penetration. Comfortable collection experience.

PATIENT SAFETY: Single-use only. Sterilised by ETO.

FDA CLEARED

ACCESSORIES

Blood Collection Needle Holders with and without safety function are compatible with our Chiravac Blood Collection Needles and our different Blood Collection Tubes. The holders are packed in PE bags of 50 or 100 pcs.

Colour Code	Size Metric	Size G/"	Article No.
Blue ●	0.60 x 38 mm	23G x 1 1/2"	CHBCN23112
Black ●	0.70 x 38 mm	22G x 1 1/2"	CHBCN22112
Green ●	0.80 x 38 mm	21G x 1 1/2"	CHBCN21112
Yellow ●	0.90 x 38 mm	20G x 1 1/2"	CHBCN20112
Pink ●	1.20 x 38 mm	18G x 1 1/2"	CHBCN18112
Blue ●	0.60 x 25 mm	23G x 1"	CHBCN23100
Black ●	0.70 x 25 mm	22G x 1"	CHBCN22100
Green ●	0.80 x 25 mm	21G x 1"	CHBCN21100
Yellow ●	0.90 x 25 mm	20G x 1"	CHBCN20100
Pink ●	1.20 x 25 mm	18G x 1"	CHBCN18100

Shelf life 5 years

Description	Pcs./Box	Pcs./Carton	Article No.
Blood Collection Needle Holder	100	1,000	CHBCNH01
Blood Collection Needle Safety Holder	50	800	CHBCNHS01

Shelf life 5 years