



Vitrex Heparin Tubes

Vitrex Heparin Tubes are coated to inhibit clotting. Heparin Tubes are used to produce plasma for biochemical testing. The Heparin Tubes are also available with Gel separator. The range of application is: clinical chemistry and trace mineral testing.

Lithium Heparin

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Plain	+ Gel
					Article No.	Article No.
LH/Lithium Heparin	1	13 x 75	Lithium Heparin	Paper	105401+ *	102701+ *
LH/Lithium Heparin	2	13 x 75	Lithium Heparin	Paper	101402+ *	102702+ *
LH/Lithium Heparin	3	13 x 75	Lithium Heparin	Paper	101403+ *	101703+ *
LH/Lithium Heparin	3	13 x 100	Lithium Heparin	Paper	101404+ **	101730+ **
LH/Lithium Heparin	4	13 x 75	Lithium Heparin	Paper	105404+ *	
LH/Lithium Heparin	4	13 x 100	Lithium Heparin	Paper	102440+ **	102704+ **
LH/Lithium Heparin	5	13 x 75	Lithium Heparin	Paper	105450+ *	
LH/Lithium Heparin	5	13 x 100	Lithium Heparin	Paper	105405+ **	102705+ **
LH/Lithium Heparin	6	13 x 100	Lithium Heparin	Paper	103406+ **	
LH/Lithium Heparin	7	13 x 100	Lithium Heparin	Paper	103407+ **	
LH/Lithium Heparin	8	16 x 100	Lithium Heparin	Paper	105408+ ***	103706+ ***
LH/Lithium Heparin	9	16 x 100	Lithium Heparin	Paper	104409+ ***	102709+ ***

Shelf life 18 months



Transit packaging

Pallet

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Sodium Heparin

Plain + Gel

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.	Article No.
NH/Sodium Heparin	1	13 x 75	Sodium Heparin	Paper	104501+ *	104601+ *
NH/Sodium Heparin	2	13 x 75	Sodium Heparin	Paper	104502+ *	104602+ *
NH/Sodium Heparin	3	13 x 75	Sodium Heparin	Paper	104530+ *	101603+ *
NH/Sodium Heparin	3	13 x 100	Sodium Heparin	Paper	104503+ **	104630+ **
NH/Sodium Heparin	4	13 x 75	Sodium Heparin	Paper	101504+ *	
NH/Sodium Heparin	4	13 x 100	Sodium Heparin	Paper	104504+ **	104604+ **
NH/Sodium Heparin	5	13 x 75	Sodium Heparin	Paper	104550+ *	
NH/Sodium Heparin	5	13 x 100	Sodium Heparin	Paper	104505+ **	102605+ **
NH/Sodium Heparin	6	13 x 100	Sodium Heparin	Paper	102506+ **	
NH/Sodium Heparin	7	13 x 100	Sodium Heparin	Paper	102507+ **	
NH/Sodium Heparin	8	16 x 100	Sodium Heparin	Paper	104509+ ***	103608+ ***
NH/Sodium Heparin	9	16 x 100	Sodium Heparin	Paper	103509+ ***	104609+ ***

Shelf life 18 months