



### Medoject® Safety Hypodermic needles

Medoject Safety Hypodermic needles have an integrated safety feature for one-handed manipulation.

| Colour Code | Size Metric  | Size G/”     | Article No. |
|-------------|--------------|--------------|-------------|
| Blue        | 0.60 x 30 mm | 23G x 1 1/4” | CHS23114    |
| Black       | 0.70 x 40 mm | 22G x 1 1/2” | CHS22112    |
| Green       | 0.80 x 40 mm | 21G x 1 1/2” | CHS21112    |
| Yellow      | 0.90 x 40 mm | 20G x 1 1/2” | CHS20112    |

Shelf life 5 years

### Medoject® Blunt fill needles

Medoject Blunt fill needles are also available with filter of 5µm for aspiration filtration procedure, which risks particle elimination during procedure.

| Colour Code | Size Metric                       | Size G/”     | Article No. |
|-------------|-----------------------------------|--------------|-------------|
| Red         | 1.20 x 40 mm Blunt fill           | 18G x 1 1/2” | CH18112BF   |
| Pink        | 1.20 x 40 mm Blunt fill w/ filter | 18G x 1 1/2” | CH18112F    |
| Pink        | 1.20 x 50 mm Blunt fill w/ filter | 18G x 2”     | CH18200F    |

Shelf life 5 years



### Medoject® Special Hypodermic needles

Medoject Special Hypodermic needles are intended for a wide special application range. There is possibility of short bevel lancet point grinding.

| Colour Code | Size Metric   | Size G/”     | Article No. |
|-------------|---------------|--------------|-------------|
| Blue        | 0.60 x 80 mm  | 23G x 3 1/8” | CH23318     |
| Green       | 0.80 x 100 mm | 21G x 4”     | CH21400     |
| Green       | 0.80 x 120 mm | 21G x 4 3/4” | CH21434     |
| Yellow      | 0.90 x 70 mm  | 20G x 2 4/5” | CH20245     |
| Pink        | 1.20 x 70 mm  | 18G x 2 4/5” | CH18245     |
| Pink        | 1.20 x 100 mm | 18G x 4”     | CH18400     |
| White       | 1.60 x 70 mm  | 16G x 2 4/5” | CH16245     |
| Blue grey   | 1.80 x 70 mm  | 15G x 2 4/5” | CH15245     |
| Pale green  | 2.00 x 70 mm  | 14G x 2 4/5” | CH14245     |
| Pale green  | 2.00 x 80 mm  | 14G x 3 1/8” | CH14318     |
| Pale green  | 2.00 x 100 mm | 14G x 4”     | CH14400     |
| Pale green  | 2.00 x 120 mm | 14G x 4 3/4” | CH14434     |

Other sizes available upon request.

Shelf life 5 years

### Inox® Ophthalmic needles

Inox Ophthalmic needles has a blunt end for ophthalmic application, which includes scratch-free edges for comfort use.

| Colour Code | Size Metric  | Size G/”   | Article No. STRAIGHT | Article No. BENT |
|-------------|--------------|------------|----------------------|------------------|
| Grey        | 0.40 x 20 mm | 27G x 4/5” | CHOP27045S           | CHOP27045B       |
| Blue        | 0.60 x 20 mm | 23G x 4/5” | CHOP23045S           | CHOP23045B       |

Shelf life 5 years