



## Blood Collection Set

The Blood Collection Sets are designed to fulfill the requirements of blood collection procedures, improve patient's comfort, help user's safety and get reliable test results. The latest products in our Blood Collection System range are of corresponding high quality and carefully selected to complete the range.

**USER VERSATILITY:** Series combination blood collection types for customer selection.

**USER CONVENIENCE:** Allows for multiple draws with a single vein puncture.

**EXCELLENT PERFORMANCE:** High quality rubber sleeve has excellent self-sealing performance.

**FLEXIBLE PACKAGING:** Single blister package, single pouch or single PE bag.

Description	Colour Code	Size	Diameter (mm)	Tube length	Pcs./ Box	Pcs./ Carton	Article No.	Article No. Safety needle
Luer Slip adapter + Intravenous needle	Green ●	21G	0.80 x 19mm	190 mm	100	2000	VBCS72221	VBCS73221S
Luer Slip adapter + Intravenous needle	Blue ●	23G	0.60 x 19mm	190 mm	100	2000	VBCS72223	VBCS73223S
Luer Slip adapter + holder + Intravenous needle	Green ●	21G	0.80 x 19mm	190 mm	50	1000	VBCS84221	VBCS84221S
Luer Slip adapter + holder + Intravenous needle	Blue ●	23G	0.60 x 19mm	190 mm	50	1000	VBCS84223	VBCS84223S

Shelf life 5 years