



## Pipette Fixed

### LOW COST

This miniature fixed volume pipette is specially designed for use with Diagnostic Kits.

### HIGH ACCURACY AND PRECISION

Each pipette is individually calibrated as per the ISO 8655 standard.

### COLOUR-CODING

The colour-coding for each pipette indicates the volume - ensuring easy identification in the laboratory environment.

### TWO STEP PLUNGER OPERATION

This unique feature makes the pipette a genuine miniature micropipette.

### PRACTICAL SIZE

At a length of 130 mm this pipette is the optimum practical size for the pipetting function.

### FULLY AUTOCLAVABLE

Autoclavable at 121°C & 15 psi.

### Technical Specifications

Colour Code	Vol. µl	Accuracy ±%	Precision CV%	Pcs.	Article No.
White ○	5	2.5	1.5	1	560013
Orange ●	10	1.2	0.8	1	561013
Black ●	20	1	0.5	1	562013
Grey ○	25	1	0.4	1	563013
Green ●	50	1	0.4	1	564013
Blue ●	100	0.8	0.3	1	565013
Red ●	200	0.8	0.3	1	566013

Shelf life 4 years