



## Allergy

For skin sensitivity testing

Allergy lancets are traditional handheld, flat and stainless steel lancet for skin sensitivity testing.

**LOW COST:** Simple, flat, steel lancet.

**HIGH QUALITY:** Hardened surgical steel.

**MINIMUM DISCOMFORT:** Modern asymmetric design.

**STERILE:** Single-use only. Sterilised by gamma irradiation.

## Steel

General Purpose Capillary Blood lancet

Steel lancets are traditional handheld, flat and stainless steel blood lancet for general purpose use.

**LOW COST:** Simple, flat, steel lancet.

**SAMPLE VOLUME:** Excellent blood flow.

**HIGH QUALITY:** Hardened surgical steel.

**MINIMUM DISCOMFORT:** Modern asymmetric design.

**STERILE:** Single-use only. Sterilised by gamma irradiation.

### Technical Specifications

Product	Description	Article No.
Steel	Flat blade sterile lancet	350213
Allergy	Flat blade sterile lancet	365213

Shelf life 5 years



**Packaging:**  
 Primary packaging: 200 pcs.  
 Transit packaging: 20,000 pcs.  
 Pallet: 720,000 pcs.