



Ultrasafe

Sterile Safety lancet

Ultrasafe lancet is intended for blood glucose and other in vitro diagnostic tests. Healthcare workers worldwide face the risk of accidental needle-stick injuries, exposing them to bloodborne illnesses including HIV and Hepatitis. Ultrasafe lancet helps reduce these risks while supporting patient comfort and precise, reliable blood collection.

Technical Specifications

Colour Code	Specification	Penetration Depth	Article No.	Indicative use
Green	17G blade	2.0 mm	393413	For capillary blood sampling
Blue	21G needle	1.8 mm	397413	Diabetes clinics/hospital wards for blood glucose monitoring
Orange	21G needle	2.2 mm	395413	Haemoglobin (blood banks), cholesterol, blood typing, coagulation, blood gas and other diagnostic tests
Pink	21G needle	2.8 mm	396413	As the Orange lancet, but for thicker adult skin
Yellow	26G needle	1.8 mm	394413	As the Orange lancet, but lower blood flow and less pain
Rose	28G needle	1.5 mm	392413	Diabetes clinics/hospital wards for blood glucose monitoring
Purple	28G needle	1.8 mm	398413	Low blood flow, less pain
Turquoise	30G needle	1.8 mm	399413	Diabetes clinics/hospital wards for blood glucose monitoring. For minimal pain

Shelf life 5 years

PRESSURE ACTIVATED

UNIQUE DESIGN:

Ergonomically optimized for better handling

SILICONE COATED NEEDLE:

Less pain and smooth penetration

EASY TARGETING OF PUNCTURE SITE



Transit packaging

Pallet

Packaging:

Primary packaging: 100 pcs.
 Transit packaging: 1,000 pcs.
 Pallet: 98,000 pcs.