



Glucose Tubes

The tubes contain a glucose stabiliser, which allows you to block the process of glycolysis for up to 24 hours at room temperature.

Range of application: Research of glucose and lactate.



Transit packaging

Pallet

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Article No.
FX/ Sodium Flouride + Potassium Oxalate	1	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	111401+ *
FX/ Sodium Flouride + Potassium Oxalate	2	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	111402+ *
FX/ Sodium Flouride + Potassium Oxalate	2	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113420+ **
FX/ Sodium Flouride + Potassium Oxalate	3	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	113430+ *
FX/ Sodium Flouride + Potassium Oxalate	3	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113403+ **
FX/ Sodium Flouride + Potassium Oxalate	4	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	112404+ *
FX/ Sodium Flouride + Potassium Oxalate	4	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113404+ **
FX/ Sodium Flouride + Potassium Oxalate	5	13 x 75	Sodium Flouride + Potassium Oxalate	Paper	113450+ *
FX/ Sodium Flouride + Potassium Oxalate	5	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113405+ **
FX/ Sodium Flouride + Potassium Oxalate	6	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113406+ **
FX/ Sodium Flouride + Potassium Oxalate	7	13 x 100	Sodium Flouride + Potassium Oxalate	Paper	113407+ **
FX/ Sodium Flouride + Potassium Oxalate	8	16 x 100	Sodium Flouride + Potassium Oxalate	Paper	113408+ ***
FX/ Sodium Flouride + Potassium Oxalate	9	16 x 100	Sodium Flouride + Potassium Oxalate	Paper	113409+ ***

Shelf life 18 months