



Heparin Tubes

Vacuum blood collection tubes with lithium heparin are used to produce plasma for biochemical testing. Gel is a special material designed to form a resistant barrier between the cellular components of blood and serum during centrifugation.

Range of application for lithium heparin: Clinical chemistry, trace mineral testing.

Range of application for lithium heparin + gel: Plasma biochemical analysis.

Lithium Heparin

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Plain	+ Gel
					Article No.	Article No.
LH/Lithium Heparin	1	13 x 75	Lithium Heparin	Paper	105401+ *	102701+ *
LH/Lithium Heparin	2	13 x 75	Lithium Heparin	Paper	101402+ *	102702+ *
LH/Lithium Heparin	3	13 x 75	Lithium Heparin	Paper	101403+ *	101703+ *
LH/Lithium Heparin	3	13 x 100	Lithium Heparin	Paper	101404+ **	101730+ **
LH/Lithium Heparin	4	13 x 75	Lithium Heparin	Paper	105404+ *	
LH/Lithium Heparin	4	13 x 100	Lithium Heparin	Paper	102440+ **	102704+ **
LH/Lithium Heparin	5	13 x 75	Lithium Heparin	Paper	105450+ *	
LH/Lithium Heparin	5	13 x 100	Lithium Heparin	Paper	105405+ **	102705+ **
LH/Lithium Heparin	6	13 x 100	Lithium Heparin	Paper	103406+ **	
LH/Lithium Heparin	7	13 x 100	Lithium Heparin	Paper	103407+ **	
LH/Lithium Heparin	8	16 x 100	Lithium Heparin	Paper	105408+ ***	103706+ ***
LH/Lithium Heparin	9	16 x 100	Lithium Heparin	Paper	104409+ ***	102709+ ***

Shelf life 18 months



Vacuum blood collection tubes with sodium heparin are used to produce plasma for biochemical testing. Gel is a special material designed to form a resistant barrier between the cellular components of blood and serum during centrifugation.

Range of application for sodium heparin: Plasma biochemical analysis.



Transit packaging

Pallet

Packaging:
 Primary packaging: 50/100 pcs.
 Transit packaging: 800/1,000 pcs.
 * Pallet: 54,000 pcs.
 ** Pallet: 48,000 pcs.
 *** Pallet: 38,400 pcs.

Sodium Heparin

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Plain		+ Gel	
					Article No.	Article No.	Article No.	Article No.
NH/Sodium Heparin	1	13 x 75	Sodium Heparin	Paper	104501+ *		104601+ *	
NH/Sodium Heparin	2	13 x 75	Sodium Heparin	Paper	104502+ *		104602+ *	
NH/Sodium Heparin	3	13 x 75	Sodium Heparin	Paper	104530+ *		101603+ *	
NH/Sodium Heparin	3	13 x 100	Sodium Heparin	Paper	104503+ **		104630+ **	
NH/Sodium Heparin	4	13 x 75	Sodium Heparin	Paper	101504+ *			
NH/Sodium Heparin	4	13 x 100	Sodium Heparin	Paper	104504+ **		104604+ **	
NH/Sodium Heparin	5	13 x 75	Sodium Heparin	Paper	104550+ *			
NH/Sodium Heparin	5	13 x 100	Sodium Heparin	Paper	104505+ **		102605+ **	
NH/Sodium Heparin	6	13 x 100	Sodium Heparin	Paper	102506+ **			
NH/Sodium Heparin	7	13 x 100	Sodium Heparin	Paper	102507+ **			
NH/Sodium Heparin	8	16 x 100	Sodium Heparin	Paper	104509+ ***		103608+ ***	
NH/Sodium Heparin	9	16 x 100	Sodium Heparin	Paper	103509+ ***		104609+ ***	

Shelf life 18 months