



## Alcohol Swabs

Alcohol Swabs are an easy and convenient way to cleanse the skin and help fight bacteria at the site of application. Use as an antiseptic skin cleanser before injections and blood sampling.

The size of the Disinfection Swab is 30x60 mm. Each Alcohol Swab is individually sealed and disposable for convenience and safety.

### Technical Specifications

Description	Size	Content	Identification No.	Pcs./Box	Article No.
Alcohol Swabs	30 x 60 mm	82% Ethanol	1U8K-YGJD-J103-YQRW	100	521413
Alcohol Swabs	30 x 60 mm	82% Ethanol + 0.5% Chlorhexidine	NQ1P-3GSC-U100-G3GE	100	522413

Shelf life 5 years

Alcohol Swabs are currently approved for sales in: Sweden, Finland, UK, Germany, France and Austria.

(Updated: 2025-10. Please consult our website: [vitrexmg.com](http://vitrexmg.com) for new updates in the list of approved countries).

Regulated according to: REGULATION (EU) No 528/2012 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 22 May 2012 concerning the making available on the market and use of biocidal products.



Transit packaging

Pallet

**Packaging:**  
 Primary packaging: 100 pcs.  
 Transit packaging: 10,000 pcs.  
 Pallet: 400,000 pcs.