



# Chirana® Insulin and Tuberculin Syringes

Chirana® Insulin and Tuberculin syringes are with integrated needle for single use. The needle cap is made of latex-free rubber gasket. The syringes are individually packed in boxes of 100 pcs. They are pyrogen-free and non-toxic.

**USER CONVENIENCE:** Needle permanently fixed to syringe body.

**USER SAFETY:** Integrated needle cap.

**COMFORT:** Grinded, precision lancet point and silicone treated tip for smooth penetration.

**PATIENT SAFETY:** Single-use only. Sterilised by ETO. Dead space is minimized to reduce waste and transmission of infectious agents. The dead space does not exceed the limits (EN ISO 8537:2016, art. 5.11.1).

**FDA CLEARED \***



Transit packaging

Pallet

**Packaging:**  
 Primary packaging: 100 pcs.  
 Secondary packaging: 2,000 pcs.  
 Pallet: 48,000 pcs.

\* FDA clearance for Insulin syringes.

## Technical Specifications

Type	Vol. ml.	Needle	Size Metric	Article No:
Tuberculin	1	29G x 1/2"	12.7 mm	CHTUB0129
Tuberculin	1	Luer		CHTUB01
Insulin U 40	1	29G x 1/2"	12.7 mm	CHINS4129
Insulin U 100	1	27G x 1/2"	12.7 mm	CHINS0127
Insulin U 100	1	29G x 5/16"	8 mm	CHINS01298
Insulin U 100	1	29G x 1/2"	12.7 mm	CHINS0129
Insulin U 100	1	30G x 1/2"	12.7 mm	CHINS0130
Insulin U 100	1	30G x 5/16"	8 mm	CHINS01308
Insulin U 100	1	31G x 4/16"	6 mm	CHINS01316
Insulin U 100	1	31G x 5/16"	8 mm	CHINS01318
Insulin U 100	1	Luer		CHINS01
Insulin U 100	0.5	29G x 1/2"	12.7 mm	CHINS00529
Insulin U 100	0.5	30G x 1/2"	12.7 mm	CHINS00530
Insulin U 100	0.5	30G x 5/16"	8 mm	CHINS005308
Insulin U 100	0.5	31G x 5/16"	8 mm	CHINS005318
Insulin U 100	0.5	31G x 4/16"	6 mm	CHINS005316
Insulin U 100	0.3	29G x 1/2"	12.7 mm	CHINS00329
Insulin U 100	0.3	30G x 1/2"	12.7 mm	CHINS00330
Insulin U 100	0.3	31G x 5/16"	8 mm	CHINS003318
Insulin U 100	0.3	31G x 4/16"	6 mm	CHINS003316

Chirana Syringes are available in many sizes.  
 For more information Contact Sales Department.

Shelf life 5 years