



Soft Blood Glucose lancets

Soft Blood Glucose lancets are intended for blood glucose measurement. The soft lancets are universal and compatible with most lancing devices especially our Compact Eject. The range includes standard soft lancets and siliconized soft lancets available in different sizes.

GOOD BLOOD FLOW

Triple ground needle point.

VERSATILE

Fit most commonly used lancing devices including the Vitrex Compact Eject lancing device.

PATIENT SAFETY:

Single-use only. Sterilised by gamma irradiation.

CONVENIENT

Boxes have dispensing feature.

NEW SILICONIZED NEEDLE

For less pain and smoother penetration.

Technical Specifications

Colour Code	Specification	Pcs./Box	Article No.
Light blue	23G (0.64 mm)	200	351213
Red	28G (0.36 mm)	200	352213
Green	30G (0.30 mm)	200	354213
Dark blue	32G (0.26 mm)	200	357213

Shelf life 5 years

Technical Specifications Siliconized Needles

Colour Code	Specification	Pcs./Box	Article No.
Red	28G (0.36 mm)	200	352413
Green	30G (0.30 mm)	200	354413

Shelf life 5 years



Transit packaging

Pallet

Packaging:
 Primary packaging: 200 pcs.
 Secondary packaging: 2,000 pcs.
 Transit packaging: 20,000 pcs.
 Pallet: 480,000 pcs.