



Accessories IV Administration

Within IV Administration for Chiraplus® Infusion set and Chirahem® Transfusion set we offer following accessories:

CHIRAPLUS® SPIKES: Intended for comfortable and safe drug manipulation.

CHIRAPLUS® INJECTION PORTS:

Offer additional drug administration when applied with other devices.

CHIRAPLUS® STOPPERS:

Give a short-term option for storage or transfer of solution.

CHIRALINE® EXTENSION LINES:

Extension tubing in 45 cm or 180 cm, where Luer lock ends are available.

CHIRAWAY® 3-WAY STOPCOCK:

Available in standard or lipid resistant material with colour-coding differentiation.

Product	Pcs./Box	Pcs./Carton	Article No.
Infusion spike Red air 0.1 µm + 5 µm solution filter	100	1000	CHIS01
Infusion spike Blue air 0.1 µm + 5 µm solution filter	100	1000	CHIS02
Infusion spike Green air 0.1 µm filter	100	1000	CHIS03
Injection port	100	5000	CHIP01
Combistopper Red	100	5000	CHCS01
Combistopper Blue	100	5000	CHCS02
Combistopper White	100	5000	CHCS03
Stopper Luer lock	100	5000	CHLL01
Chiraline® 450	120	960	CHL318045
Chiraline® 1800	85	680	CHL318180
Chiraway® Standard	50	1000	CHW001
Chiraway® Lipid resistant	50	1000	CHW002

Shelf life 5 years