



## Chiraflex® IV Cannulas and Chiraflex® Safety IV Cannulas

Chiraflex® IV Cannulas and Chiraflex® Safety IV Cannulas are intended for peripheral vascular application. The catheter is made from Teflon material and available in Safety. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. They are packed in boxes of 50 pcs, cartons of 1000 pcs.

### USER CONVENIENCE:

Ompatible with all standard medical devices for wide application.

**USER SAFETY:** Safety cannula is available.

**PATIENT SAFETY:** Single-use only. Sterilised by ETO.

| Colour Code | Size G/" | Size Metric | Article No. With Port With Wings | Article No. Without Port With Wings | Article No. Without Port Without Wings | Article No. SAFETY With Port With Wings |
|-------------|----------|-------------|----------------------------------|-------------------------------------|--|---|
| Purple ●    | 26G      | 0.60 x 19mm | CHFX0102619                      | CHFX0202619                         | CHFX0302619                            | CHFXS0102619                            |
| Yellow ●    | 24G      | 0.70 x 19mm | CHFX0102419                      | CHFX0202419                         | CHFX0302419                            | CHFXS0102419                            |
| Blue ●      | 22G      | 0.90 x 25mm | CHFX0102225                      | CHFX0202225                         | CHFX0302225                            | CHFXS0102225                            |
| Pink ●      | 20G      | 1.10 x 32mm | CHFX0102032                      | CHFX0202032                         | CHFX0302032                            | CHFXS0102032                            |
| Green ●     | 18G      | 1.30 x 45mm | CHFX0101845                      | CHFX0201845                         | CHFX0301845                            | CHFXS0101845                            |
| White ○     | 17G      | 1.50 x 45mm | CHFX0101745                      | CHFX0201745                         | CHFX0301745                            | CHFXS0101745                            |
| Grey ●      | 16G      | 1.70 x 45mm | CHFX0101645                      | CHFX0201645                         | CHFX0301645                            | CHFXS0101645                            |
| Orange ●    | 14G      | 2.10 x 45mm | CHFX0101445                      | CHFX0201445                         | CHFX0301445                            | CHFXS0101445                            |

Shelf life 5 years