



## Chiraflex® Mandrins

Chiraflex® Mandrins are intended for longer interruptions in infusion and has a screw onto IV cannula. The mandrins are compatible with Chiraflex® IV Cannulas and Chiraflex® Safety IV Cannulas. They are colour-coded as per Chiraflex® IV Cannula colour-coding.

Colour Code	Size	Pcs./Box	Pcs./Carton	Article No.
Blue	22G	50	1000	CHFX20022
Pink	20G	50	1000	CHFX20020
Green	18G	50	1000	CHFX20018
White	17G	50	1000	CHFX20017
Grey	16G	50	1000	CHFX20016
Orange	14G	50	1000	CHFX20014

Shelf life 5 years

## Chiraflex® Scalp vein sets

Chiraflex® Scalp vein sets are intended for shorter time infusion. The set includes a stainless steel needle with a length of 20 mm and a 30 cm plastic tube.

Colour Code	Size	Diameter (mm)	Pcs./Box	Pcs./Carton	Article No.
Grey	27G x 3/4"	0.40	100	2000	CH260001127
Orange	25G x 3/4"	0.50	100	2000	CH260001125
Blue	23G x 3/4"	0.60	100	2000	CH260001123
Black	22G x 3/4"	0.70	100	2000	CH260001122
Green	21G x 3/4"	0.80	100	2000	CH260001121
Yellow	20G x 3/4"	0.90	100	2000	CH260001120
Ivory	19G x 3/4"	1.10	100	2000	CH260001119
Pink	18G x 3/4"	1.20	100	2000	CH260001118

Shelf life 5 years